

AURELIA® Quest 2.2®

Examination Gloves



92935	(XS)	ø 200	10 x 200
92936	(S)	ø 200	10 x 200
92937	(M)	ø 200	10 x 200
92938	(L)	ø 200	10 x 200
92939	(XL)	ø 200	10 x 200



Nitrile



Powder Free



White



3.4+/- 0.2 g



Finger-Textured




5 Years

Premium White Nitrile examination gloves that offer increased tensile strength and protection over Latex and Vinyl.

Nitrile examination gloves offer increased comfort and flexibility.

- 2.2mil Powder-Free Nitrile
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Latex Free
- AQL 1.5

Sizes & Dimensions



XS	70 ± 10 mm	min. 240 mm
S	80 ± 10 mm	min. 240 mm
M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	≥ 110 mm	min. 240 mm

Single wall thickness



- ① min. 0.11 mm
- ② min. 0.05 mm
- ③ min. 0.05 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 14.0 MPa	min. 14.0 MPa
Elongation	min. 500 %	min. 400 %
Strength at Break	min. 6N	min. 6N

Quality Signs & Certificates

- 89/686/EEC, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 2-3, EN 388
- EN 1186: Food Contact Approved (SGS)

