



AURELIA® Robust Plus®

Examination Gloves



63885	(XS)	ø 100	10 x 100
63886	(S)	ø 100	10 x 100
63887	(M)	ø 100	10 x 100
63888	(L)	ø 100	10 x 100
63889	(XL)	ø 100	10 x 100



Nitrile



Powder Free



Regular Blue



5.8 +/- 0.2 g



Micro-Textured




5 Years

Premium Quality Blue Nitrile examination gloves that offer greater tensile strength and protection than both Latex and Vinyl.

Nitrile examination gloves offer extra comfort and flexibility and longer length cuff for increased protection.

- 5.0mil Powder-Free Nitrile
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Latex Free
- AQL 1.5

Sizes & Dimensions



Size	Width (mm)	Length (mm)
XS	70 ± 10 mm	min. 290 mm
S	80 ± 10 mm	min. 290 mm
M	95 ± 10 mm	min. 290 mm
L	110 ± 10 mm	min. 290 mm
XL	≥ 110 mm	min. 290 mm

Single wall thickness



- ① min. 0.13 mm
- ② min. 0.10 mm
- ③ min. 0.07 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 14.0 MPa	min. 14.0 MPa
Elongation	min. 500 %	min. 400 %
Strength at Break	min. 6N	min. 6N

Quality Signs & Certificates

- 89/686/EEC, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 2-3, EN 388
- EN 1186: Food Contact Approved (SGS)

