

AURELIA® Robust®

Examination Gloves



| | | | |
|--------------|------|-------|----------|
| 93895 | (XS) | ø 100 | 10 x 100 |
| 93896 | (S) | ø 100 | 10 x 100 |
| 93897 | (M) | ø 100 | 10 x 100 |
| 93898 | (L) | ø 100 | 10 x 100 |
| 93899 | (XL) | ø 100 | 10 x 100 |



Nitrile



Powder Free



Regular Blue



5.0 +/- 0.2 g



Micro-Textured




5 Years

Premium Quality Blue Nitrile examination gloves that offer greater tensile strength and protection than both Latex and Vinyl.

Nitrile examination gloves offer extra comfort and flexibility.

- 4.5mil Powder-Free Nitrile
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Latex Free
- AQL 1.5

Sizes & Dimensions



| | | |
|-----------|-------------|-------------|
| XS | 70 ± 10 mm | min. 240 mm |
| S | 80 ± 10 mm | min. 240 mm |
| M | 95 ± 10 mm | min. 240 mm |
| L | 110 ± 10 mm | min. 240 mm |
| XL | ≥ 110 mm | min. 240 mm |

Single wall thickness



- ① min. 0.12 mm
- ② min. 0.09 mm
- ③ min. 0.05 mm

Physical properties

| | Unaged | Aged |
|-------------------|---------------|---------------|
| Tensile Strength | min. 14.0 MPa | min. 14.0 MPa |
| Elongation | min. 500 % | min. 400 % |
| Strength at Break | min. 6N | min. 6N |

Quality Signs & Certificates

- 89/686/EEC, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 2-3, EN 388
- EN 1186: Food Contact Approved (SGS)

