









AURELIA® Sonic 200®

Examination Gloves



9377A5 (XS)	ø 200	10 x 200
9377A6 (S)	ø 200	10 x 200
9377A7 (M)	ø 200	10 x 200
9377A8 (L)	ø 200	10 x 200
9377A9 (XL)	ø 200	10 x 200


-  **Nitrile**
-  **Powder Free**
-  **Cobalt Blue**
-  **Finger-Textured**
-  **3.2+/- 0.2 g**
-  **5 Years**

Quality Cobalt Blue Nitrile examination gloves that offer increased tensile strength and protection over Latex and Vinyl.

Incredible Ultra-Thin Nitrile examination gloves specifically designed to offer increased comfort, flexibility & sensitivity.

- 2.2mil Powder-Free Nitrile
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Latex Free
- AQL 1.5

Sizes & Dimensions



XS	70 ± 10 mm	min. 240 mm
S	80 ± 10 mm	min. 240 mm
M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	≥ 110 mm	min. 240 mm

Single wall thickness



- ① min. 0.11 mm
- ② min. 0.05 mm
- ③ min. 0.05 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 14.0 MPa	min. 14.0 MPa
Elongation	min. 500 %	min. 400 %
Strength at Break	min. 6N	min. 6N

Quality Signs & Certificates

- 89/686/EEC, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 2-3, EN 388
- EN 1186: Food Contact Approved (SGS)

