

AURELIA® Vibrant®

Examination Gloves



98225	(XS)	ø 100	10 x 100
98226	(S)	ø 100	10 x 100
98227	(M)	ø 100	10 x 100
98228	(L)	ø 100	10 x 100
98229	(XL)	ø 100	10 x 100



Latex



Powder Free



Natural



5.7 +/- 0.2 g



Micro Textured




3 Years

Premium Natural Latex examination gloves that offer superior strength, comfort, dexterity and sensitivity.

- Chlorinated Powder-Free Latex
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Micro Textured
- AQL 1.5

Sizes & Dimensions



Size	Width (mm)	Length (mm)
XS	70 ± 10 mm	min. 240 mm
S	80 ± 10 mm	min. 240 mm
M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	≥ 110 mm	min. 240 mm

Single wall thickness



- ① min. 0.13 mm
- ② min. 0.10 mm
- ③ min. 0.07 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 18.0 MPa	min. 14.0 MPa
Elongation	min. 650 %	min. 500 %
Strength at Break	min. 6N	min. 6N

Quality Signs & Certificates

- 89/686/EEC, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 2-3, EN 388
- EN 1186: Food Contact Approved (SGS)

