



AURELIA® ROBUST PLUS

4.5mil, Extended Cuff, Nitrile Powder Free Examination Gloves

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk
+441733 615215

Supermax Healthcare Ltd,
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN

aureliagloves.co.uk



10 x 100




XS	S	M	L	XL
63885	63886	63887	63889	63899

Sizes & Dimensions

Single Wall Thickness

Physical Properties



XS	70 - 79 mm	290 ± 10 mm
S	80 - 89 mm	290 ± 10 mm
M	90 - 99 mm	290 ± 10 mm
L	100 - 109 mm	290 ± 10 mm
XL	110 - 119 mm	290 ± 10 mm

- ① min. 0.13 mm
- ② min. 0.10 mm
- ③ min. 0.07 mm

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N



Nitrile



Powder Free



5 Years



Micro



5.8g +/- 0.2g



Regular Blue

