



# AURELIA® UNIQUE TPE CLEAR

Powder Free TPE Examination Gloves

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk

+441733 615215

Supermax Healthcare Ltd,  
Unit 12-16 Titan Drive  
Fengate, Peterborough  
PE1 5XN

[aureliagloves.co.uk](http://aureliagloves.co.uk)






Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 455 1-4 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.

 10 x 200

**S M L XL**

49226 49227 49228 49229

Sizes & Dimensions		Single wall thickness	Physical Properties													
																
S	103+2 mm	min. 270 mm	① min. 0.06 mm													
M	105+2 mm	min. 270 mm	② min. 0.06 mm													
L	115+2 mm	min. 270 mm	③ min. 0.05 mm													
XL	> 120 mm	min. 270 mm														
				<table border="1"> <thead> <tr> <th></th> <th>Unaged</th> <th>Aged</th> </tr> </thead> <tbody> <tr> <td>Tensile Strength</td> <td>Min. 14.0 MPa</td> <td>Min. 14.0 MPa</td> </tr> <tr> <td>Elongation</td> <td>Min. 400%</td> <td>Min. 400%</td> </tr> <tr> <td>Strength at Break</td> <td>Min. 3.0N</td> <td>Min. 3.0N</td> </tr> </tbody> </table>		Unaged	Aged	Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa	Elongation	Min. 400%	Min. 400%	Strength at Break	Min. 3.0N	Min. 3.0N
	Unaged	Aged														
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa														
Elongation	Min. 400%	Min. 400%														
Strength at Break	Min. 3.0N	Min. 3.0N														



TPE



Powder Free



3 Years



Smooth



2.2g +/- 0.2g



Clear

