



# AURELIA® VIBRANT

## Premium Natural Latex Disposable Gloves

### AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk  
+441733 615215

Supermax Healthcare Ltd,  
Unit 12-16 Titan Drive  
Fengate, Peterborough  
PE1 5XN

[aureliagloves.co.uk](http://aureliagloves.co.uk)

<b>Regulations</b>	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
<b>Quality Standards</b>	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
<b>Harmonised Standards</b>	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.



10 x 100



<b>XS</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>
98265	98226	98227	98228	98229

#### Sizes & Dimensions

#### Single wall thickness

#### Physical Properties

XS 70+10 mm	min. 240 mm
S 80+10 mm	min. 240 mm
M 95+10 mm	min. 240 mm
L 110+10 mm	min. 240 mm
XL > 110 mm	min. 240 mm



- ① min. 0.13 mm
- ② min. 0.10 mm
- ③ min. 0.08 mm

	Unaged	Aged
Tensile Strength	Min. 18.0 MPa	Min. 14.0 MPa
Elongation	Min. 650%	Min. 500%
Strength at Break	Min. 6N	Min. 6N



Latex



Powder Free



3 Years



Micro



5.7g +/- 0.2g



Natural

