



AURELIA® VELOCITY

Premium Natural Latex Disposable Gloves

AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk
+441733 615215

Supermax Healthcare Ltd,
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN

aureliagloves.co.uk

Regulations	PPE Regulation (EU) 2016/42 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.



10 x 100




XS	S	M	L	XL
88265	88266	88267	88268	88269

Sizes & Dimensions

Single wall thickness

Physical Properties



XS	70+10 mm	min. 240 mm
S	80+10 mm	min. 240 mm
M	95+10 mm	min. 240 mm
L	110+10 mm	min. 240 mm
XL	> 110 mm	min. 240 mm

①	min. 0.15 mm
②	min. 0.12 mm
③	min. 0.08 mm

	Unaged	Aged
Tensile Strength	Min. 18.0 MPa	Min. 14.0 MPa
Elongation	Min. 650%	Min. 500%
Strength at Break	Min. 6N	Min. 6N



Latex



Powder Free



3 Years



Micro



6.5g +/- 0.2g



Natural

