

# **AURELIA® BLUSH**

2.5mil, Nitrile Powder Free Examination Gloves

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC	
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems	
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.	

## AURELIA® GLOVES A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk +441733 615215

Supermax Healthcare Ltd, Unit 12-16 Titan Drive Fengate, Peterborough PE1 5XN

aureliagloves.co.uk



















78885 78886 78887 78888 78889

#### **Sizes & Dimensions**

#### Single wall thickness

### **Physical Properties**



XS 70+10 mm 80+10 mm 95+10 mm

L 110+10 mm XL > 110 mm



min. 240 mm min. 240 mm min. 240 mm min. 240 mm min. 240 mm



1 min. 0.10 mm 2 min. 0.05 mm

3 min. 0.05 mm

	ogou	1.900
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N



**Nitrile** 



**Powder Free** 





**Finger Textured** 



3.4g +/- 0.2g





