



AURELIA® DELIGHT CLEAR PF

Premium Powder Free Vinyl Examination Gloves

AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk

+441733 615215

Supermax Healthcare Ltd,
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN

aureliagloves.co.uk



Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 455 1-4 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.



10 x 100

S

M

L

XL

38226

38227


38228

38229

Sizes & Dimensions

Single wall thickness

Physical Properties

		
S	80+10 mm	min. 240 mm
M	95+10 mm	min. 240 mm
L	110+10 mm	min. 240 mm
XL	> 110 mm	min. 240 mm



- ① min. 0.07 mm
- ② min. 0.05 mm
- ③ min. 0.03 mm

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.6N	Min. 3.6N



Vinyl



Powder Free



3 Years



Smooth



4.5g +/- 0.2g



Clear

