



# AURELIA® PERFORM

2.5mil, Nitrile Powder Free Examination Gloves



<b>Regulations</b>	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
<b>Quality Standards</b>	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
<b>Harmonised Standards</b>	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1184: 1997. Materials and articles in contact with foodstuffs.

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk

+441733 615215

Supermax Healthcare Ltd,  
Unit 12-16 Titan Drive  
Fengate, Peterborough  
PE1 5XN

[aureliagloves.co.uk](http://aureliagloves.co.uk)



10 x 100

**XS S M L XL**

92395 92396 92397 92398 92399

### Sizes & Dimensions



XS 70+10 mm  
S 80+10 mm  
M 95+10 mm  
L 110+10 mm  
XL > 110 mm



min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm

### Single wall thickness



① min. 0.10 mm  
② min. 0.05 mm  
③ min. 0.05 mm

### Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N



Nitrile



Powder Free



5 Years



Finger  
Textured



3.4g +/- 0.2g



Teal Green

