

AURELIA® ROBUST PLUS

4.5mil, Nitrile Powder Free Examination Gloves

AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk
+441733 615215

Supermax Healthcare Ltd,
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN

aureliagloves.co.uk



| | |
|----------------------|--|
| Regulations | PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC |
| Quality Standards | ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1184: 1997. Materials and articles in contact with foodstuffs. |



10 x 100


XS S M L XL

63885 63886 63887 63889 63899

Sizes & Dimensions

Single wall thickness

Physical Properties



XS 70+10 mm
S 80+10 mm
M 95+10 mm
L 110+10 mm
XL > 110 mm



min. 290 mm
min. 290 mm
min. 290 mm
min. 290 mm
min. 290 mm



① min. 0.13 mm
② min. 0.10 mm
③ min. 0.07 mm

| | Unaged | Aged |
|-------------------|---------------|---------------|
| Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa |
| Elongation | Min. 500% | Min. 400% |
| Strength at Break | Min. 6N | Min. 6N |



Nitrile



Powder Free



5 Years



Textured



5.8g +/- 0.2g



Sky Blue

