



AURELIA® DELIGHT BLUE PD

Blue Vinyl Powdered Examination Gloves

AURELIA® GLOVES

sales@supermax.co.uk
01733 615215

Supermax Healthcare Ltd
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN



Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186:2002 Materials and articles in contact with foodstuffs.




10 x 100

S	M	L	XL
38896	38897	38898	38899

Sizes & Dimensions

Single Wall Thickness

Physical Properties



S	80±10 mm
M	95±10 mm
L	110±10 mm
XL	> 110 mm



min. 240 mm
min. 240 mm
min. 240 mm
min. 240 mm



- ① min. 0.05mm
- ② min. 0.05mm
- ③ min. 0.04mm

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.6N	Min. 3.6N



Vinyl



Powdered



3 Years



Smooth



4.5g +/- 0.2g



Blue

