



# AURELIA® DELIGHT BLUE PF

Blue Vinyl Powder Free Examination Gloves

## AURELIA® GLOVES

sales@supermax.co.uk  
01733 615215




Supermax Healthcare Ltd  
Unit 12-16 Titan Drive  
Fengate, Peterborough  
PE1 5XN




Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186:2002 Materials and articles in contact with foodstuffs.

 10 x 100

S	M	L	XL
38996	38997	38998	38999

Sizes & Dimensions		Single Wall Thickness	Physical Properties	
			Unaged	Aged
S 80±10 mm	min. 240 mm	① min. 0.06mm	Tensile Strength	Min. 11.0 MPa
M 95±10 mm	min. 240 mm	② min. 0.06mm	Elongation	Min. 350%
L 110±10 mm	min. 240 mm	③ min. 0.05mm	Strength at Break	Min. 3.6N
XL > 110 mm	min. 240 mm			

-  Vinyl
-  Powder Free
-  3 Years
-  Smooth
-  4.5g +/- 0.2g
-  Blue

