



AURELIA® DELIGHT CLEAR PF AURELIA® GLOVES

Vinyl Powder Free Examination Gloves

sales@supermax.co.uk
01733 615215

Supermax Healthcare Ltd
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN



Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods BS EN 455 1-4 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.

10 x 100

S	M	L	XL
38226	38227	38228	38229

Sizes & Dimensions	Single Wall Thickness	Physical Properties													
 S 80±10 mm M 95±10 mm L 110±10 mm XL > 110 mm	 min. 240 mm min. 240 mm min. 240 mm min. 240 mm	 ① min. 0.05mm ② min. 0.05mm ③ min. 0.04mm	<table border="1"> <thead> <tr> <th></th> <th>Unaged</th> <th>Aged</th> </tr> </thead> <tbody> <tr> <td>Tensile Strength</td> <td>Min. 11.0 MPa</td> <td>Min. 11.0 MPa</td> </tr> <tr> <td>Elongation</td> <td>Min. 350%</td> <td>Min. 350%</td> </tr> <tr> <td>Strength at Break</td> <td>Min. 3.6N</td> <td>Min. 3.6N</td> </tr> </tbody> </table>		Unaged	Aged	Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa	Elongation	Min. 350%	Min. 350%	Strength at Break	Min. 3.6N	Min. 3.6N
	Unaged	Aged													
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa													
Elongation	Min. 350%	Min. 350%													
Strength at Break	Min. 3.6N	Min. 3.6N													



Vinyl



Powder Free



3 Years



Smooth



4.5g +/- 0.2g



Clear

