



AURELIA® DELIGHT LITE

Vinyl Powder Free Examination Gloves

AURELIA® GLOVES

sales@supermax.co.uk
01733 615215

Supermax Healthcare Ltd
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN



Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems. ISO 13485: 2016 Medical Devices Quality Management Systems.
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 455 1&3 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.



10 x 100

S

M

L

XL

32226

32227


32228

32229

Sizes & Dimensions

Single Wall Thickness

Physical Properties



S 85±5 mm
M 95±5 mm
L 105±5 mm
XL > 115 mm



min. 240 mm
min. 240 mm
min. 240 mm
min. 240 mm



① min. 0.05mm
② min. 0.05mm
③ min. 0.04mm

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.6N	Min. 3.6N



Vinyl



Powder Free



5 Years



Smooth



4.0g +/- 0.2g



Clear

