



AURELIA® IGNITE

7mil, Nitrile Powder Free Examination Gloves

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems. ISO 13485: 2016 Medical Devices Quality Management Systems.
Harmonised Standards	BS EN ISO 374 - 1 Protective gloves against dangerous chemicals, Type B JKPT. BS EN ISO 374 - 5 Protective gloves against microorganisms (Virus). BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186:2002 Materials and articles in contact with foodstuffs.

AURELIA® GLOVES

sales@supermax.co.uk
 01733 615215
 Supermax Healthcare Ltd
 Unit 12-16 Titan Drive
 Fengate, Peterborough
 PE1 5XN

10 x 100
 10 x 90*

M	L	*XL	*XXL
97887	97888	97889	97880

Sizes & Dimensions	Single Wall Thickness	Physical Properties													
<p>M 95±10 mm</p> <p>L 110±10 mm</p> <p>XL 120±10 mm</p> <p>XXL > 120 mm</p>	<p>min. 240 mm</p> <p>min. 240 mm</p> <p>min. 240 mm</p> <p>min. 240 mm</p> <p>min. 240 mm</p>	<p>① min. 0.25 mm</p> <p>② min. 0.20 mm</p> <p>③ min. 0.09 mm</p>	<table border="1"> <thead> <tr> <th></th> <th>Unaged</th> <th>Aged</th> </tr> </thead> <tbody> <tr> <td>Tensile Strength</td> <td>Min. 14.0 MPa</td> <td>Min. 14.0 MPa</td> </tr> <tr> <td>Elongation</td> <td>Min. 500%</td> <td>Min. 400%</td> </tr> <tr> <td>Strength at Break</td> <td>Min. 6N</td> <td>Min. 6N</td> </tr> </tbody> </table>		Unaged	Aged	Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa	Elongation	Min. 500%	Min. 400%	Strength at Break	Min. 6N	Min. 6N
	Unaged	Aged													
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa													
Elongation	Min. 500%	Min. 400%													
Strength at Break	Min. 6N	Min. 6N													



Nitrile



Powder Free



5 Years



Diamond Textured



8.4g +/- 0.2g



Orange

