



AURELIA® UNIQUE TPE BLUE

Blue TPE Powder Free Examination Gloves

AURELIA® GLOVES

sales@supermax.co.uk
01733 615215

Supermax Healthcare Ltd
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN



Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 455 1 & 2 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.

10 x 200

S	M	L	XL
48226	48227	48228	48229

Sizes & Dimensions

Single Wall Thickness

Physical Properties

S 103±2 mm	min. 270 mm
M 105±2 mm	min. 270 mm
L 115±2 mm	min. 270 mm
XL > 120 mm	min. 270 mm

① min. 0.06mm
② min. 0.06mm
③ min. 0.05mm

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 400%	Min. 400%
Strength at Break	Min. 3.0N	Min. 3.0N



TPE



Powder Free



3 Years



Smooth



2.2g +/- 0.2g



Blue

