



# AURELIA® SONIC 100

2.2mil, Nitrile Powder Free Examination Gloves

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk

+441733 615215

Supermax Healthcare Ltd,

Unit 12-16 Titan Drive


Fengate, Peterborough

PE1 5XN

[aureliagloves.co.uk](http://aureliagloves.co.uk)



|                             |  |
|-----------------------------|--|
| <b>Regulations</b>          | PPE Regulation (EU) 2016/425<br>Medical Device Regulation (EU) 2017/745, 93/42/EEC   |
| <b>Quality Standards</b>    | ISO 9001:2015 Quality Management Systems<br>ISO 13485: 2016 Medical Devices Quality Management Systems   |
| <b>Harmonised Standards</b> | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms.<br>BS EN 455 1-4 Medical gloves for single use.<br>EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods.<br>BS EN 1184: 1997. Materials and articles in contact with foodstuffs. |



| Sizes & Dimensions  | Single wall thickness  | Physical Properties  |  |        |      |                  |               |               |            |           |           |                   |         |         |
|---|--|--|--|--------|------|------------------|---------------|---------------|------------|-----------|-----------|-------------------|---------|---------|
| <p>S 80+10 mm<br/>M 95+10 mm<br/>L 110+10 mm<br/>XL &gt; 110 mm</p> | <p>min. 240 mm<br/>min. 240 mm<br/>min. 240 mm<br/>min. 240 mm</p> | <p>① min. 0.13 mm<br/>② min. 0.09 mm<br/>③ min. 0.05 mm</p>  |  |        |      |                  |               |               |            |           |           |                   |         |         |
|   |  | <table border="1"> <thead> <tr> <th></th> <th>Unaged</th> <th>Aged</th> </tr> </thead> <tbody> <tr> <td>Tensile Strength</td> <td>Min. 14.0 MPa</td> <td>Min. 14.0 MPa</td> </tr> <tr> <td>Elongation</td> <td>Min. 500%</td> <td>Min. 400%</td> </tr> <tr> <td>Strength at Break</td> <td>Min. 6N</td> <td>Min. 6N</td> </tr> </tbody> </table> |  | Unaged | Aged | Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa | Elongation | Min. 500% | Min. 400% | Strength at Break | Min. 6N | Min. 6N |
|   | Unaged   | Aged   |  |        |      |                  |               |               |            |           |           |                   |         |         |
| Tensile Strength  | Min. 14.0 MPa  | Min. 14.0 MPa  |  |        |      |                  |               |               |            |           |           |                   |         |         |
| Elongation  | Min. 500%  | Min. 400%  |  |        |      |                  |               |               |            |           |           |                   |         |         |
| Strength at Break   | Min. 6N  | Min. 6N  |  |        |      |                  |               |               |            |           |           |                   |         |         |



Nitrile



Powder Free



5 Years



Finger Textured



3.0g +/- 0.2g



Cobalt Blue

