



AURELIA® ROBUST 9.0

Heavy Duty 9 Newton Finger Textured Blue Nitrile Disposable Gloves

AURELIA® ROBUST 9.0

Heavy Duty 9 Newton Finger Textured Blue Nitrile Disposable Gloves

- Nitrile**
- Powder Free**
- 5.5g**
- 5.0 mil**
- Finger Textured**
- Screen Compatible**

Applications

10 x 100

S	M	L	XL
96896	96897	96898	96899

MD 0598

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 9N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
240 x 125 x 85 mm

Outer Case Dimensions:
430 x 255 x 248 mm

EAN 13 Barcodes:
Inner
S: 955-500210-2876
M: 955-500210-2883
L: 955-500210-2890
XL: 955-500210-2906

Outer
S: 955-500210-2873
M: 955-500210-2880
L: 955-500210-2897
XL: 955-500210-2903


Features and Benefits

Premium Blue Nitrile: The strongest blue nitrile glove in our range, it carries the industry known, and popular Robust name.

9 Newton Tensile Strength: Designed to provide enhanced durability, puncture resistance, and comfort, offering superior strength over the industry standard 6 newton glove.

Finger-Textured: Featuring a textured finish on the fingers only for improved grip, dexterity and tactile sensitivity.

Single Wall Thickness Sizes & Dimensions



- ① min. 0.13 mm
- ② min. 0.10 mm
- ③ min. 0.08 mm



Size	Length	Min. Thickness
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie

