





AURELIA® ABSOLUTE 100


Finger Textured Black Nitrile Disposable Gloves - Class I Medical Device

**Nitrile**

**Powder Free**

**3.6g**

**3.2 mil**

**Finger Textured**

**Screen Compatible**

 **10 x 100**

XS	S	M	L	XL
9899A5	9899A6	9899A7	9899A8	9899A9

**MD**

**CE**
0598


**AQL 1.5**

**VIRUS**

**JKT**

**LATEX**

Applications



AURELIA® ABSOLUTE 100

Finger Textured Black Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
235 x 125 x 60 mm

Outer Case Dimensions:
315 x 260 x 245 mm

EAN 13 Barcodes:
Inner

XS: 955-500211-9775
S: 955-500211-9782
M: 955-500211-9799
L: 955-500211-9805
XL: 955-500211-9812

Outer

XS: 1955-500211-9772
S: 1955-500211-9789
M: 1955-500211-9796
L: 1955-500211-9802
XL: 1955-500211-9819


Features and Benefits

All Rounder: Versatile and reliable, suitable for a wide range of tasks and environments, offering balanced protection and comfort.

Dexterity: Enhanced precision and control for intricate tasks, ensuring optimal performance in various sectors.

Finger-Textured: Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

Single Wall Thickness Sizes & Dimensions



① min. 0.10 mm
② min. 0.06 mm
③ min. 0.05 mm



XS	70 - 79 mm	min. 240
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie




@aureliaglovesuk