






AURELIA® AMAZING  
Finger Textured Ice Blue Nitrile Disposable Gloves - Class I Medical Device




Nitrile




Powder Free




3.0g



2.2 mil



Finger Textured



Screen Compatible





10 x 300




\*10 x 270


| XS    | S     | M     | L     | *XL   |
|-------|-------|-------|-------|-------|
| 92885 | 92886 | 92887 | 92888 | 92889 |




0598




AQL 1.5



BS EN ISO 374-5:2016  
VIRUS



BS EN ISO 374-1:2016  
JKT



LATEX

Applications

AURELIA® AMAZING  
Finger Textured Ice Blue Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

|                   | Unaged        | Aged          |
|-------------------|---------------|---------------|
| Tensile Strength  | Min. 14.0 MPa | Min. 14.0 MPa |
| Elongation        | Min. 500%     | Min. 400%     |
| Strength at Break | Min. 6N       | Min. 6N       |

Regulations

|                      |   |
|----------------------|---|
| Regulations          | PPE Regulation (EU) 2016/425<br>Medical Device Regulation (EU) 2017/745, 93/42/EEC  |
| Quality Standards    | ISO 9001:2015 Quality Management Systems<br>ISO 13485: 2016 Medical Devices Quality Management Systems  |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms<br>EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods.<br>BS EN 1186: 2002. Materials and articles in contact with foodstuffs.<br>BS EN 455 1-4 Medical gloves for single use. |

Inner Box Dimensions:  
240 x 125 x 85 mm

Outer Case Dimensions:  
440 x 255 x 248 mm

EAN 13 Barcodes:  
Inner  
XS: 955-500210-2807  
S: 955-500210-2814  
M: 955-500210-2821  
L: 955-500210-2838  
XL: 955-500210-2845

Outer  
XS: 1955-500210-2804  
S: 1955-500210-2811  
M: 1955-500210-2828  
L: 1955-500210-2835  
XL: 1955-500210-2842


Features and Benefits

**Lightweight 3g:** Designed for comfort and dexterity without compromising on strength, these gloves are ideal for tasks requiring precision and agility.


**Finger-Textured:** Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

**300 Gloves Per Box:** Providing an increased quantity per box for optimal efficiency and prolonged use, perfect for environments requiring frequent glove changes.

Single Wall Thickness      Sizes & Dimensions



- ① min. 0.09 mm
- ② min. 0.05 mm
- ③ min. 0.05 mm



|    |              |          |
|----|--------------|----------|
| XS | 70 - 79 mm   | min. 240 |
| S  | 80 - 89 mm   | min. 240 |
| M  | 90 - 99 mm   | min. 240 |
| L  | 100 - 109 mm | min. 240 |
| XL | 110 - 119 mm | min. 240 |



**UK & Europe HQ**  
Supermax Healthcare Limited  
Units 12-16 Titan Drive,  
Fengate, Peterborough,  
PE1 5XN, UK

+441733 615215  
sales@supermax.co.uk  
www.aureliagloves.co.uk



**Europe**  
Supermax Healthcare (Europe) Limited  
Units 1001-1012  
Building 1000, Gateway Business Park,  
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724  
sales@supermax.ie  
www.aureliagloves.ie



  
@aureliaglovesuk