



AURELIA® BLUSH


Finger Textured Pink Nitrile Disposable Gloves - Class I Medical Device




Nitrile




Powder Free




3.4g




2.5 mil



Finger Textured





Screen Compatible




10 x 200


XS	S	M	L	XL
78885	78886	78887	78888	78889




0598






AQL 1.5



BS EN ISO 374-5:2016



BS EN ISO 374-1:2016 Type C



Applications



AURELIA® BLUSH

Finger Textured Pink Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
240 x 125 x 85 mm

Outer Case Dimensions:
430 x 255 x 248 mm

EAN 13 Barcodes:
Inner
XS: 955-500210-5242
S: 955-500210-5259
M: 955-500210-5266
L: 955-500210-5273
XL: 955-500210-5280

Outer
XS: 1955-500210-5249
S: 1955-500210-5256
M: 1955-500210-5263
L: 1955-500210-5270
XL: 1955-500210-5287


Features and Benefits

Bright Pink: A popular and in-demand colour, offering high visibility and a vibrant look while maintaining exceptional performance and protection.


Finger-Textured: Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

200 Gloves Per Box: Packaged to maximise storage facilities, cost-effectiveness and efficiency, ensuring ample supply for busy high use settings.

Single Wall Thickness Sizes & Dimensions



① min. 0.10 mm
② min. 0.05 mm
③ min. 0.05 mm



XS	70 - 79 mm	min. 240
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie





@aureliaglovesuk