



AURELIA® BOLD MAX

Heavy Duty Diamond Textured Black Nitrile Disposable Gloves

Nitrile

Powder Free

7.5g

6.0 mil

Diamond Textured

Screen Compatible

10 x 50

M	L	XL	XXL
9789A7	9789A8	9789A9	9789A0

MD

CE

AQL 1.5

VIRUS

JKPT

LATEX

0598

Applications

AURELIA® BOLD MAX

Heavy Duty Diamond Textured Black Nitrile Disposable Gloves

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 9N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
240 x 130 x 70 mm

Outer Case Dimensions:
360 x 270 x 255 mm

EAN 13 Barcodes:
Inner

M: 955-500211-3063
L: 955-500211-3070
XL: 955-500211-3087
XXL: 955-500211-3094

Outer

M: 1955-500211-3060
L: 1955-500211-3077
XL: 1955-500211-3084
XXL: 1955-500211-3091

Features and Benefits

50 Gloves Per Box: Our heavy duty cost-efficient offering, designed for convenient storage and packing.

Heavy Duty 7.5g: The second heaviest nitrile disposable glove in our range that thrives in demanding and tough environments.

Diamond Textured: Engineered with a diamond textured grip for enhanced gripping and control, providing excellent performance in wet and oily conditions.

Single Wall Thickness Sizes & Dimensions

① min. 0.20 mm
② min. 0.17 mm
③ min. 0.09 mm

M	90 - 99 mm	min. 240 mm
L	100 - 109 mm	min. 240 mm
XL	110 - 119 mm	min. 240 mm
XXL	120 - 129 mm	min. 240 mm



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork,T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie



@aureliaglovesuk