



AURELIA® BOLD

Heavy Duty Micro Textured Black Nitrile Disposable Gloves

Nitrile

Powder Free

5.5g

5.0 mil

Micro Textured

Screen Compatible

10 x 100

S	M	L	XL	XXL
73996	73997	73998	73999	73990

MD

0598

AQL 1.5

BS EN ISO 374-5: 2016 VIRUS

BS EN ISO 374-1: 2016 JKPT

Applications

AURELIA® BOLD

Heavy Duty Micro Textured Black Nitrile Disposable Gloves

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
240 x 125 x 85 mm

Outer Case Dimensions:
430 x 255 x 248 mm

EAN 13 Barcodes:
Inner
S: 955-500210-5204
M: 955-500210-5211
L: 955-500210-5228
XL: 955-500210-5235
XXL: 955-500210-5341

Outer
S: 1955-500210-5201
M: 1955-500210-5218
L: 1955-500210-5225
XL: 1955-500210-5232
XXL: 1955-500210-5348

Features and Benefits

Premium Black Nitrile: The demand for black nitrile gloves is at an all-time high, driven by increasing consumer and industry popularity.

Rigorous Performance: Designed to withstand demanding conditions, offering superior durability and protection for intense applications.

Micro-Textured: Designed for enhanced grip and control, providing excellent tactile sensitivity and precision for a variety of tasks across multiple industries.

Single Wall Thickness Sizes & Dimensions

min. 0.13 mm
min. 0.11 mm
min. 0.08 mm

S	80 - 89 mm	min. 240 mm
M	90 - 99 mm	min. 240 mm
L	100 - 109 mm	min. 240 mm
XL	110 - 119 mm	min. 240 mm
XXL	120 - 129 mm	min. 240 mm

UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk

Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie



@aureliaglovesuk