



AURELIA® DELIGHT CLEAR PF

Powder Free Clear Vinyl Disposable Gloves




Vinyl



Powder Free



4.5g



Smooth Textured



Screen Compatible



10 x 100

S	M	L	XL
38226	38227	38228	38229



0598



AQL 1.5



BS EN ISO 374-5: 2016



BS EN ISO 374-1: 2016/Type C



Applications





AURELIA® DELIGHT CLEAR PF

Powder Free Clear Vinyl Disposable Gloves

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.5N	Min. 3.5N

Regulations

Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Innecase Dimensions:
235 x 125 x 60 mm

Outercase Dimensions:
325 x 255 x 250 mm

EAN 13 Barcodes:
Inner
S: 955-500210-6072
M: 955-500210-6089
L: 955-500210-6096
XL: 955-500210-6102

Outer
S: 5056730301271
M: 5056730301288
L: 5056730301295
XL: 5056730301301


Features and Benefits

Powder Free Eliminates the risk of powder-related irritation and contamination, providing a clean and comfortable user experience.

Food Contact Approved Meets safety standards for direct food contact, ensuring hygienic handling and preventing contamination in food preparation environments.


Smooth-Textured Ensures a comfortable fit and effortless handling, offering sleek surface for easy grip and reduced friction.

Single Wall Thickness



- 1 min. 0.07 mm
- 2 min. 0.05 mm
- 3 min. 0.03 mm

Sizes & Dimensions



- S 80+10 mm
- M 95+10 mm
- L 110+10 mm
- XL >110 mm



- min. 240 mm
- min. 240 mm
- min. 240 mm
- min. 240 mm



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork,T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie

