






AURELIA® IGNITE
Heavy Duty Diamond Textured Orange Nitrile Disposable Gloves




Nitrile




Powder Free




8.4g



7.0 mil



Diamond Textured



Screen Compatible





10 x 100




10 x 90


M	L	*XL	*XXL
97887	97888	97889	97880




0598




AQL 1.5




BS EN ISO 374-5: 2016



BS EN ISO 374-1: 2016/Type B







JKPT






LATEX

Applications









AURELIA® IGNITE
Heavy Duty Diamond Textured Orange Nitrile Disposable Gloves

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 9N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
240 x 130 x 95 mm

Outer Case Dimensions:
490 x 270 x 255 mm

EAN 13 Barcodes:
Inner
M: 9 55500 2110208
L: 9 55500 2110215
XL: 9 55500 2110222
XXL: 9 55500 2110239

Outer
M: 19 55500 211 0205
L: 19 55500 211 0212
XL: 19 55500 211 0229
XXL: 19 55500 211 0236

Features and Benefits

High-Vis Orange: Featuring a bright, high-visibility colour for enhanced safety and visibility in low-light conditions, making it ideal for a range of industrial and automotive environments.

Heavy Duty 8.4g: The strongest and heaviest weighted disposable glove in our range, designed to withstand the toughest working conditions.

Diamond Textured: Engineered with a diamond textured grip for enhanced gripping and control, providing excellent performance in wet and oily conditions.

Single Wall Thickness Sizes & Dimensions



- 1 min. 0.25 mm
- 2 min. 0.20 mm
- 3 min. 0.09 mm



M	90 - 99 mm	min. 240 mm
L	100 - 109 mm	min. 240 mm
XL	110 - 119 mm	min. 240 mm
XXL	120 - 129 mm	min. 240 mm



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie

