



AURELIA® PERFORM

Finger Textured Teal Green Nitrile Disposable Gloves - Class I Medical Device

Nitrile

Powder Free

3.4g

2.5 mil

Finger Textured

Screen Compatible

10 x 200

XS	S	M	L	XL
92395	92396	92397	92398	92399

MD

CE

AQL 1.5

BS EN ISO 374-5:2016

BS EN ISO 374-1:2016/Type C

VIRUS

LATEX

Applications



AURELIA® PERFORM

Finger Textured Teal Green Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
240 x 125 x 85 mm

Outer Case Dimensions:
430 x 255 x 248 mm

EAN 13 Barcodes:
Inner
XS: 955-500210-2746
S: 955-500210-2753
M: 955-500210-2760
L: 955-500210-2777
XL: 955-500210-2784

Outer
XS: 1955-500210-2743
S: 1955-500210-2750
M: 1955-500210-2767
L: 1955-500210-2774
XL: 1955-500210-2781


Features and Benefits

Teal Green: Popular in the medical and dental markets for its distinctive colour, which aids in easy identification and enhances visibility during procedures.

Finger-Textured: Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

200 Gloves per Box: Packaged to maximise storage facilities, cost-effectiveness and efficiency, ensuring ample supply for busy high use settings.

Single Wall Thickness Sizes & Dimensions



- 1 min. 0.10 mm
- 2 min. 0.05 mm
- 3 min. 0.05 mm



XS	70 - 79 mm	min. 240
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240




UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie



@aureliaglovesuk