





AURELIA® QUEST 2.2


Finger Textured White Nitrile Disposable Gloves - Class I Medical Device




Nitrile




Powder Free




3.0g



2.2 mil



Finger Textured



Screen Compatible

10 x 200

| S | M | L | XL |
|-------|-------|-------|-------|
| 92936 | 92937 | 92938 | 92939 |




0598



AQL 1.5




BS EN ISO 374-5: 2016



BS EN ISO 374-1: 2016 Type B



VIRUS



JKT



LATEX

Applications



AURELIA® QUEST 2.2

Finger Textured White Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

| | Unaged | Aged |
|-------------------|---------------|---------------|
| Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa |
| Elongation | Min. 500% | Min. 400% |
| Strength at Break | Min. 6N | Min. 6N |

Regulations

| | |
|----------------------|--|
| Regulations | PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC |
| Quality Standards | ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs. |

Inner Box Dimensions:
235 x 125 x 70 mm

Outer Case Dimensions:
365 x 260 x 244 mm


EAN 13 Barcodes:
Inner
XS: 955-500210-5143
S: 955-500210-5150
M: 955-500210-5167
L: 955-500210-5174
XL: 955-500210-5181

Outer
XS: 955-500210-5140
S: 955-500210-5157
M: 955-500210-5164
L: 955-500210-5171
XL: 955-500210-5188


Features and Benefits

- Lightweight 3g:** Designed for comfort and dexterity without compromising on strength, these gloves are ideal for tasks requiring precision and agility.
- Finger-Textured:** Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.
- 200 Gloves per Box:** Packaged to maximise storage facilities, cost-effectiveness and efficiency, ensuring ample supply for busy high use settings.

Single Wall Thickness Sizes & Dimensions



① min. 0.09 mm
② min. 0.05 mm
③ min. 0.05 mm



| | | |
|----|--------------|----------|
| XS | 70 - 79 mm | min. 240 |
| S | 80 - 89 mm | min. 240 |
| M | 90 - 99 mm | min. 240 |
| L | 100 - 109 mm | min. 240 |
| XL | 110 - 119 mm | min. 240 |




UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie





@aureliaglovesuk