



AURELIA® QUEST

Finger Textured White Nitrile Disposable Gloves - Class I Medical Device

Nitrile

Powder Free

3.8g

3.5 mil

Finger Textured

Screen Compatible

10 x 100

XS	S	M	L	XL
92895	92896	92897	92898	92899

MD

CE

AQL 1.5

BS EN ISO 374-5:2016

BS EN ISO 374-1:2016 Type B

VIRUS

KPT

LATEX

Applications



AURELIA® QUEST

Finger Textured White Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:

245 x 130 x 60 mm

Outer Case Dimensions:

315 x 270 x 255 mm

EAN 13 Barcodes:

Inner

XS: 955-500210-2425
S: 955-500210-2432
M: 955-500210-2449
L: 955-500210-2456
XL: 955-500210-2463


Outer

XS: 1955-500210-2422
S: 1955-500210-2439
M: 1955-500210-2446
L: 1955-500210-2453
XL: 1955-500210-2460

Features and Benefits

- All Rounder:** Versatile and reliable, suitable for a wide range of tasks and environments, offering balanced protection and comfort.
- Instrument Handling:** Provides a secure grip and excellent tactile sensitivity, ensuring precision and control with tools.
- Finger-Textured:** Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

Single Wall Thickness Sizes & Dimensions



- ① min. 0.12 mm
- ② min. 0.06 mm
- ③ min. 0.05 mm



XS	70 - 79 mm	min. 240
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie

