





AURELIA® REFRESH


Honeycomb Textured Polymer Coated Mint Scented Latex Disposable Gloves




Latex




Powder Free




7.2g



Honeycomb



Screen Compatible



Mint Scented

10 x 100


XS	S	M	L	XL
99225	99226	99227	99228	99229



0598




AQL 1.5



VIRUS



Applications



AURELIA® REFRESH

Honeycomb Textured Polymer Coated Mint Scented Latex Disposable Gloves

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 18.0 MPa	Min. 14.0 MPa
Elongation	Min. 650%	Min. 500%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Innerbox Dimensions:
240 x 125 x 85 mm

Outerbox Dimensions:
430 x 255 x 248 mm

EAN 13 Barcodes:
Inner

XS: 955-500210-1558
S: 955-500210-1565
M: 955-500210-1572
L: 955-500210-1589
XL: 955-500210-1596


Outer

XS: 1955-500210-1555
S: 1955-500210-1562
M: 1955-500210-1579
L: 1955-500210-1586
XL: 1955-500210-1593

Features and Benefits

- Mint-Scented:** Provides a mint scented ingredient that enhances comfort and freshness, reducing odor and improving the overall user experience.
- Teal Green:** Popular in the medical and dental markets for its distinctive colour, which aids in easy identification and enhances visibility during procedures.
- Honeycomb Textured:** Engineered with a unique diamond texture for superior grip and control, providing excellent performance in wet and oily conditions.

Single Wall Thickness



- 1 min. 0.15 mm
- 2 min. 0.12 mm
- 3 min. 0.06 mm

Sizes & Dimensions


	Length	Width
XS	70 - 79 mm	min. 240 mm
S	80 - 89 mm	min. 240 mm
M	90 - 99 mm	min. 240 mm
L	100 - 109 mm	min. 240 mm
XL	110 - 119 mm	min. 240 mm


UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk


Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie




@aureliaglovesuk