





# AURELIA® ROBUST 9.0


Heavy Duty 9 Newton Finger Textured Blue Nitrile Disposable Gloves




Nitrile




Powder Free




5.5g




5.0 mil



Finger Textured





Screen Compatible




10 x 100


S	M	L	XL
96896	96897	96898	96899




0598




AQL 1.5



VIRUS




JKPT



LATEX

## Applications



# AURELIA® ROBUST 9.0

Heavy Duty 9 Newton Finger Textured Blue Nitrile Disposable Gloves

## Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 9N	Min. 6N

## Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:  
240 x 125 x 85 mm

Outer Case Dimensions:  
430 x 255 x 248 mm

EAN 13 Barcodes:  
Inner  
S: 955-500210-5365  
M: 955-500210-5372  
L: 955-500210-5389  
XL: 955-500210-5396

Outer  
S: 1955-500210-2873  
M: 1955-500210-2880  
L: 1955-500210-2897  
XL: 1955-500210-2903


## Features and Benefits

**Premium Blue Nitrile:** The strongest blue nitrile glove in our range, it carries the industry known, and popular Robust name.

**9 Newton Tensile Strength:** Designed to provide enhanced durability, puncture resistance, and comfort, offering superior strength over the industry standard 6 newton glove.

**Finger-Textured:** Featuring a textured finish on the fingers only for improved grip, dexterity and tactile sensitivity.

## Single Wall Thickness Sizes & Dimensions



- 1 min. 0.13 mm
- 2 min. 0.10 mm
- 3 min. 0.08 mm



S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240



**UK & Europe HQ**  
Supermax Healthcare Limited  
Units 12-16 Titan Drive,  
Fengate, Peterborough,  
PE1 5XN, UK

+441733 615215  
sales@supermax.co.uk  
www.aureliagloves.co.uk



**Europe**  
Supermax Healthcare (Europe) Limited  
Units 1001-1012  
Building 1000, Gateway Business Park,  
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724  
sales@supermax.ie  
www.aureliagloves.ie

