





AURELIA® ROBUST PLUS


Micro Textured Blue Nitrile Disposable Gloves with an Extended Cuff




Nitrile




Powder Free




5.7g




4.5 mil



Micro Textured





Screen Compatible




10 x 100


S	M	L	XL
63886	63887	63888	63889




0598




AQL 1.5




BS EN ISO 374-5:2016




BS EN ISO 374-1:2016/Type B



VIRUS




KPT



LATEX

Applications



AURELIA® ROBUST PLUS

Micro Textured Blue Nitrile Disposable Gloves with an Extended Cuff

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
270 x 125 x 85 mm

Outer Case Dimensions:
440 x 255 x 280 mm

EAN 13 Barcodes:

Inner

XS: 955-500210-5358
S: 955-500210-5365
M: 955-500210-5372
L: 955-500210-5389
XL: 955-500210-5396

Outer

XS: 1955-500210-5355
S: 1955-500210-5362
M: 1955-500210-5379
L: 1955-500210-5386
XL: 1955-500210-5393


Features and Benefits

Extended Cuff: Designed with an extended cuff (290mm) for added wrist protection, ensuring a secure fit and enhanced coverage.

Micro-Textured: Designed for enhanced grip and control, providing excellent tactile sensitivity and precision for a variety of tasks across multiple industries.

Added Splash Protection: The additional extension of the cuff ensures safety and compliance in environments where spills are a concern.

Single Wall Thickness Sizes & Dimensions



- 1 min. 0.13 mm
- 2 min. 0.10 mm
- 3 min. 0.07 mm



S	80 - 89 mm	290 ± 10 mm
M	90 - 99 mm	290 ± 10 mm
L	100 - 109 mm	290 ± 10 mm
XL	110 - 119 mm	290 ± 10 mm



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie

