






AURELIA® ROBUST
Heavy Duty Micro Textured Blue Nitrile Disposable Gloves




Nitrile




Powder Free




5.0g



4.5 mil





Micro Textured




Screen Compatible

10 x 100


| XS | S | M | L | XL | XXL |
|-------|-------|-------|-------|-------|-------|
| 93895 | 93896 | 93897 | 93898 | 93899 | 93890 |




0598




AQL 1.5



VIRUS



KPT



LATEX

Applications



AURELIA® ROBUST
Heavy Duty Micro Textured Blue Nitrile Disposable Gloves

Physical Properties

| | Unaged | Aged |
|-------------------|---------------|---------------|
| Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa |
| Elongation | Min. 500% | Min. 400% |
| Strength at Break | Min. 6N | Min. 6N |

Regulations

| | |
|----------------------|--|
| Regulations | PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC |
| Quality Standards | ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs. |

Inner Box Dimensions:
240 x 125 x 85 mm

Outer Case Dimensions:
430 x 255 x 248 mm

EAN 13 Barcodes:
Inner
XS: 955-500210-2203
S: 955-500210-2210
M: 955-500210-2227
L: 955-500210-2234
XL: 955-500210-2241
XXL: 955-500210-2982

Outer
XS: 955-500210-2200
S: 955-500210-2217
M: 955-500210-2224
L: 955-500210-2231
XL: 955-500210-2248
XXL: 955-500210-2989


Features and Benefits

Industry Favourite: Trusted for reliability and performance, widely known and used across various sectors.


Full Size Range: Available in a variety of sizes to ensure a proper fit for all users, enhancing comfort and effectiveness.

Micro-Textured: Designed for enhanced grip and control, providing excellent tactile sensitivity and precision for a variety of tasks across multiple industries.


Single Wall Thickness Sizes & Dimensions



1 min. 0.12 mm
2 min. 0.09 mm
3 min. 0.05 mm




| | | |
|-----|--------------|----------|
| XS | 70 - 79 mm | min. 240 |
| S | 80 - 89 mm | min. 240 |
| M | 90 - 99 mm | min. 240 |
| L | 100 - 109 mm | min. 240 |
| XL | 110 - 119 mm | min. 240 |
| XXL | 120 - 129 mm | min. 240 |


UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk


Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie




@aureliaglovesuk