



# AURELIA® SONIC 100

Finger Textured Cobalt Blue Nitrile Disposable Gloves - Class I Medical Device

Nitrile

Powder Free

3.0g

2.2 mil

Finger Textured

Screen Compatible

10 x 100

XS	S	M	L	XL
9377B5	9377B6	9377B7	9377B8	9377B9

MD

0598

AQL 1.5

BS EN ISO 374-5:2016

BS EN ISO 374-1:2016/Type B

VIRUS

JKT

LATEX

### Applications

# AURELIA® SONIC 100

Finger Textured Cobalt Blue Nitrile Disposable Gloves - Class I Medical Device

### Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

### Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:  
200 x 110 x 60 mm

Outer Case Dimensions:  
310 x 225 x 210 mm

EAN 13 Barcodes:  
Inner  
XS: 9555002107055  
S: 9555002107062  
M: 9 55500 210 7079  
L: 9555002107086  
XL: 9555002107093

Outer  
XS: 19 55500 210 7076  
S: 19 55500 210 7069  
M: 19 55500 210 7076  
L: 19 55500 210 7083  
XL: 19 55500 210 7090


### Features and Benefits

**Lightweight 3g:** Designed for comfort and dexterity without compromising on strength, these gloves are ideal for tasks requiring precision and agility.


**Finger-Textured:** Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

**Increased Dexterity:** Enhanced precision and control for intricate tasks, ensuring optimal performance in a variety of sectors.

### Single Wall Thickness      Sizes & Dimensions



- 1 min. 0.09 mm
- 2 min. 0.05 mm
- 3 min. 0.05 mm



XS	70 - 79 mm	min. 240
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240



**UK & Europe HQ**  
Supermax Healthcare Limited  
Units 12-16 Titan Drive,  
Fengate, Peterborough,  
PE1 5XN, UK

+441733 615215  
sales@supermax.co.uk  
www.aureliagloves.co.uk



**Europe**  
Supermax Healthcare (Europe) Limited  
Units 1001-1012  
Building 1000, Gateway Business Park,  
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724  
sales@supermax.ie  
www.aureliagloves.ie



@aureliaglovesuk