




AURELIA® SONIC 200  
Finger Textured Cobalt Blue Nitrile Disposable Gloves - Class I Medical Device




Nitrile




Powder Free




3.0g



2.2 mil





Finger Textured




Screen Compatible

10 x 200


XS	S	M	L	XL
9377A5	9377A6	9377A7	9377A8	9377A9




0598





AQL 1.5



BS EN ISO 374-5:2016 VIRUS



BS EN ISO 374-1:2016 JKT



LATEX

Applications



AURELIA® SONIC 200  
Finger Textured Cobalt Blue Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:  
235 x 125 x 70 mm

Outer Case Dimensions:  
360 x 244 x 260 mm


EAN 13 Barcodes:  
Inner  
XS: 9 55500 210 566 2  
S: 9 55500 210 567 9  
M: 9 55500 210 568 6  
L: 9 55500 210 569 3  
XL: 9 55500 210 570 9

Outer  
XS: 19 55500 210 566 9  
S: 19 55500 210 567 6  
M: 19 55500 210 568 3  
L: 19 55500 210 568 3  
XL: 19 55500 210 570 6


Features and Benefits

- Latex Free:** Contains no natural latex rubber; making it perfect for a variety of medical environments.
- Screen Compatible:** Products are screen compatible with tablets and mobile devices - if contaminated; please clean screens after use.
- 200 Gloves per Box:** Packaged to maximise storage facilities, cost-effectiveness and efficiency, ensuring ample supply for busy high use settings.

Single Wall Thickness      Sizes & Dimensions



① min. 0.09 mm  
② min. 0.05 mm  
③ min. 0.05 mm



XS	70 - 79 mm	min. 240
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240



**UK & Europe HQ**  
Supermax Healthcare Limited  
Units 12-16 Titan Drive,  
Fengate, Peterborough,  
PE1 5XN, UK

+441733 615215  
sales@supermax.co.uk  
www.aureliagloves.co.uk



**Europe**  
Supermax Healthcare (Europe) Limited  
Units 1001-1012  
Building 1000, Gateway Business Park,  
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724  
sales@supermax.ie  
www.aureliagloves.ie





@aureliaglovesuk