



AURELIA® TRANSFORM 100
Finger Textured Ice Blue Nitrile Disposable Gloves - Class I Medical Device

Nitrile

PF

3.6g

3.2 mil

Finger Textured

Screen Compatible

10 x 100

XS S M L XL
9889A5 9889A6 9889A7 9889A8 9889A9

MD

0598

AQL 1.5

BS EN ISO 374-5: 2016

BS EN ISO 374-1: 2016/Type B

Applications

AURELIA® TRANSFORM 100
Finger Textured Ice Blue Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
235 x 125 x 60 mm

Outer Case Dimensions:
315 x 260 x 245 mm

EAN 13 Barcodes:
Inner
XS: 955-500210-2920
S: 955-500210-2937
M: 955-500210-2944
L: 955-500210-2951
XL: 955-500210-2968

Outer
XS: 955-500210-2927
S: 955-500210-2934
M: 955-500210-2941
L: 955-500210-2958
XL: 955-500210-2965

Features and Benefits

- All Rounder:** Versatile and reliable, suitable for a wide range of tasks and environments, offering balanced protection and comfort.
- Instrument Handling:** Provides a secure grip and excellent tactile sensitivity, ensuring precision and control with tools.
- Finger-Textured:** Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

Single Wall Thickness Sizes & Dimensions

1 min. 0.10 mm
2 min. 0.06 mm
3 min. 0.05 mm

XS	70 - 79 mm	min. 240
S	80 - 89 mm	min. 240
M	90 - 99 mm	min. 240
L	100 - 109 mm	min. 240
XL	110 - 119 mm	min. 240

UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk

Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie



@aureliaglovesuk