






AURELIA® VIBRANT
Micro Textured Chlorinated Natural Latex Disposable Gloves




Latex




Powder Free




5.7g



5.0 mil





Micro Textured




Screen Compatible

10 x 100


XS	S	M	L	XL
98225	98226	98227	98228	98229




0598







AQL 1.5




BS EN ISO 374-5: 2016



BS EN ISO 374-1: 2016/Type C



Applications



AURELIA® VIBRANT
Heavy Duty Micro Textured Black Nitrile Disposable Gloves

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 18.0 MPa	Min. 14.0 MPa
Elongation	Min. 650%	Min. 500%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Outer case Dimensions:
365 x 260 x 244 mm

Innerrcase Dimensions:
235 x 125 x 70 mm

EAN 13 Barcodes:
Inner
XS: 955-500210-1053
S: 955-500210-1060
M: 955-500210-1077
L: 955-500210-1084
XL: 955-500210-1091

Outer
XS: 1955-500210-1050
S/7: 1955-500210-1067
M/8: 1955-500210-1074
L/9: 1955-500210-1081
XL/10: 1955-500210-1098


Features and Benefits

Chlorinated: Enhanced for a smoother, easier fit and improved resistance to tears, offering added durability and comfort.

Natural Latex: Offers exceptional elasticity and comfort, derived from natural rubber for a snug fit and reliable performance.

Micro-Textured: Designed for enhanced grip and control, providing excellent tactile sensitivity and precision for a variety of tasks across multiple industries.

Single Wall Thickness



- 1 min. 0.13 mm
- 2 min. 0.10 mm
- 3 min. 0.07 mm

Sizes & Dimensions


	XS	S	M	L	XL
Hand Length	70 - 79 mm	80 - 89 mm	90 - 99 mm	100 - 109 mm	110 - 119 mm
Hand Width	min. 240 mm	min. 240 mm	min. 240 mm	min. 240 mm	min. 240 mm

UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk

Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie



@aureliaglovesuk