



AURELIA® ZERO

Finger Textured Accelerator Free Ice Blue Nitrile Disposable Gloves - Class I Medical Device

Nitrile

Powder Free

3.6g

3.2 mil

Finger Textured

Screen Compatible

10 x 200

XS	S	M	L	XL
92775	92776	92777	92778	92779

MD

CE

AQL 1.5


BS EN ISO 374-5:2016

BS EN ISO 374-1:2016/Type C

VIRUS

LATEX

Applications



AURELIA® ZERO

Finger Textured Accelerator Free Ice Blue Nitrile Disposable Gloves - Class I Medical Device

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Regulations

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Inner Box Dimensions:
240 x 125 x 85 mm

Outer Case Dimensions:
430 x 255 x 248 mm

EAN 13 Barcodes:
Inner Box
XS: 955-500210-5617
S: 955-500210-5624
M: 955-500210-5631
L: 955-500210-5648
XL: 955-500210-5655

Outer
XS: 955-500210-5614
S: 955-500210-5621
M: 955-500210-5638
L: 955-500210-5645
XL: 955-500210-5652

Features and Benefits

Accelerator Free: Minimises the risk of allergic reactions and skin irritation, enhancing safety and comfort for users.

200 Gloves per Box: Packaged to maximise storage facilities, cost-effectiveness and efficiency, ensuring ample supply for busy high use settings.

Finger-Textured: Featuring a textured finish on the fingers only for improved grip, dexterity, tactile sensitivity.

Single Wall Thickness Sizes & Dimensions

XS	70-79 mm	min. 240
S	80-89 mm	min. 240
M	90-99 mm	min. 240
L	100-109 mm	min. 240
XL	110-119 mm	min. 240

1	min. 0.10 mm
2	min. 0.06 mm
3	min. 0.05 mm




UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork,T23 F2H7,Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie



@aureliaglovesuk