



# AURELIA® VINTAGE

Powdered Latex Disposable Gloves

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk

+441733 615215

Supermax Healthcare Ltd,  
Unit 12-16 Titan Drive  
Fengate, Peterborough  
PE1 5XN

[aureliagloves.co.uk](http://aureliagloves.co.uk)

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.




10 x 100

**XS S M L XL**  
 28825 28826 28827 28828 28829

### Sizes & Dimensions

### Single wall thickness

### Physical Properties



XS 70+10 mm  
S 80+10 mm  
M 95+10 mm  
L 110+10 mm  
XL > 110 mm



min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm



① min. 0.08 mm  
② min. 0.08 mm  
③ min. 0.06 mm

	Unaged	Aged
Tensile Strength	Min. 18.0 MPa	Min. 14.0 MPa
Elongation	Min. 650%	Min. 500%
Strength at Break	Min. 6N	Min. 6N



Latex



Powdered



3 Years



Smooth



5.5g +/- 0.2g



Natural

