

# **AURELIA® ROBUST 9.0**

5mil, 9 Newton, Nitrile Powder Free Examination Gloves

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

## AURELIA® GLOVES

sales@supermax.co.uk +441733 615215

Supermax Healthcare Ltd, Unit 12-16 Titan Drive Fengate, Peterborough PE1 5XN

aureliagloves.co.uk

96895 96896 96897 96898 96899











min. 240 mm












#### Sizes & Dimensions

### Single wall thickness

- ① min. 0.13 mm
- 2 min. 0.10 mm
- 3 min. 0.08 mm

#### **Physical Properties**





XS 70+10 mm

L 110+10 mm

XL > 110 mm

80+10 mm

95+10 mm













