



AURELIA® DELIGHT CLEAR PD

Powdered Clear Vinyl Disposable Gloves




Vinyl



Powdered



4.5g



Smooth Textured



Screen Compatible

10 x 100

S	M	L	XL
38826	38827	38828	38829



BS EN ISO 374-5: 2016



VIRUS



KPT



LATEX

Applications





AURELIA® DELIGHT CLEAR PD

Powdered Clear Vinyl Disposable Gloves

Physical Properties

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.5N	Min. 3.5N

Regulations

Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Innecase Dimensions:
235 x 125 x 60 mm

Outercase Dimensions:
325 x 255 x 250 mm


EAN 13 Barcodes:
Inner
S: 955-500210-6027
M: 955-500210-6034
L: 955-500210-6041
XL: 955-500210-6058

Outer
S: 195-55002-106024
M: 195-55002-106031
L: 195-55002-106048
XL: 195-55002-106109

Features and Benefits

- Food Contact Approved:** Meets safety standards for direct food contact, ensuring hygienic handling and preventing contamination in food preparation environments.
- Slightly Powdered:** Provides a light dusting to facilitate easier donning and removal, while minimising moisture and stickiness.
- Smooth-Textured:** Ensures a comfortable fit and effortless handling, offering a sleek surface for easy grip and reduced friction.

Single Wall Thickness



- ① min. 0.07 mm
- ② min. 0.05 mm
- ③ min. 0.03 mm

Sizes & Dimenensions

	80+10 mm	min. 240 mm
S	80+10 mm	min. 240 mm
M	95+10 mm	min. 240 mm
L	110+10 mm	min. 240 mm
XL	>110 mm	min. 240 mm



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie





@aureliaglovesuk