



AURELIA® BOLD

5mil, Nitrile Powder Free Examination Gloves

AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk
+441733 615215

Supermax Healthcare Ltd,
Unit 12-16 Titan Drive
Fengate, Peterborough
PE1 5XN

aureliagloves.co.uk

| | |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Regulations | PPE Reugulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC |
| Quality Standards | ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs. |

 10 x 100














S M L XL XXL
 73996 73997 73998 73999 73990

| Sizes & Dimensions | Single wall thickness | Physical Properties | | | | | | | | | | | | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------|------|------------------|---------------|---------------|------------|-----------|-----------|-------------------|---------|---------|
|  <p>XS 70+10 mm S 80+10 mm M 95+10 mm L 110+10 mm XL > 110 mm</p> |  <p>min. 240 mm min. 240 mm min. 240 mm min. 240 mm min. 240 mm</p> |  <p>① min. 0.13 mm ② min. 0.10 mm ③ min. 0.08 mm</p> <table border="1"> <thead> <tr> <th></th> <th>Unaged</th> <th>Aged</th> </tr> </thead> <tbody> <tr> <td>Tensile Strength</td> <td>Min. 14.0 MPa</td> <td>Min. 14.0 MPa</td> </tr> <tr> <td>Elongation</td> <td>Min. 500%</td> <td>Min. 400%</td> </tr> <tr> <td>Strength at Break</td> <td>Min. 6N</td> <td>Min. 6N</td> </tr> </tbody> </table> | | Unaged | Aged | Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa | Elongation | Min. 500% | Min. 400% | Strength at Break | Min. 6N | Min. 6N |
| | Unaged | Aged | | | | | | | | | | | | |
| Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa | | | | | | | | | | | | |
| Elongation | Min. 500% | Min. 400% | | | | | | | | | | | | |
| Strength at Break | Min. 6N | Min. 6N | | | | | | | | | | | | |



Nitrile



Powder Free



5 Years



Micro



5.5g +/- 0.2g



Black

