



# AURELIA® DELIGHT BLUE PF

Premium Powder Free Vinyl Examination Gloves

**AURELIA® GLOVES**  
A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk  
+441733 615215

Supermax Healthcare Ltd,  
Unit 12-16 Titan Drive  
Fengate, Peterborough  
PE1 5XN




[aureliagloves.co.uk](http://aureliagloves.co.uk)



Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186:2002 Materials and articles in contact with foodstuffs.

 10 x 100

S	M	L	XL
38996	38997	38998	38999

	Sizes & Dimensions		Single wall thickness	Physical Properties	
				Unaged	Aged
					
S	80+10 mm	min. 240 mm	① min. 0.06 mm	Min. 11.0 MPa	Min. 11.0 MPa
M	95+10 mm	min. 240 mm	② min. 0.06 mm		
L	110+10 mm	min. 240 mm	③ min. 0.05 mm	Min. 350%	Min. 350%
XL	> 110 mm	min. 240 mm		Min. 3.6N	Min. 3.6N



Vinyl



Powder Free



3 Years



Smooth



4.5g +/- 0.2g



Blue

