



# AURELIA® ROBUST 9.0

5mil, 9 Newton, Nitrile Powder Free Examination Gloves

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.co.uk

+441733 615215

Supermax Healthcare Ltd,

Unit 12-16 Titan Drive

Fengate, Peterborough

PE1 5XN

[aureliagloves.co.uk](http://aureliagloves.co.uk)



10 x 100

XS S M L XL


96895 96896 96897 96898 96899



### Sizes & Dimensions

### Single wall thickness

### Physical Properties



XS 70+10 mm  
S 80+10 mm  
M 95+10 mm  
L 110+10 mm  
XL > 110 mm



min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm



① min. 0.13 mm  
② min. 0.10 mm  
③ min. 0.08 mm

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 9N	Min. 9N



Nitrile



Powder Free



5 Years



Finger Textured



5.5g +/- 0.2g



Ice Blue

